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The RRBS-seq experiment will be conducted by DNA Methylation Profiling (RRBS-seq Service reduced-representation bisulfite sequencing) in a total of 21 octopus samples.

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DNA concentration of the samples will be measured using the Qubit® dsDNA BR Assay Kit (Thermo Fisher Scientific). DNA quality of the samples will be assessed with the Fragment Analyzer™ and the DNF-488 High Sensitivity genomic DNA Analysis Kit (Agilent).

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RRBS libraries will be prepared using the Premium Reduced Representation Bisulfite Sequencing (RRBS) Kit (Diagenode Cat# C02030033). 100ng of genomic DNA were used to start library preparation for each sample. Following library preparation, samples were pooled together by 10. PCR clean-up after the final library amplification was performed using a 1.45x beads:sample ratio of Agencourt® AMPure® XP (Beckman Coulter).

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DNA concentration of the pools will be measured using the Qubit® dsDNA HS Assay Kit (Thermo Fisher Scientific). The profile of the pools will be checked using the DNF-474 NGS fragment kit on a Fragment Analyzer (Agilent).

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RRBS library pools were sequenced on a NovaSeq6000 (Illumina) using 50 bp paired-end read sequencing (PE50).

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Customized RRBS bioinformatic analysis following Klughammer et al., 20215 for analysis without reference genome.

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Comparison of methylation status of CpG nucleotides between sample groups.

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Gene ontology, enrichment and pathway analysis.

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